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Please read carefully.

EXPLANATION OF SYMBOLS

The following signs and symbols on items, packaging and in the instruction manual bear important information:

 This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

 **WARNING**
These warning notes must be observed to prevent any injury to the user.

 **CAUTION**
These notes must be observed to prevent any damage to the device.

 **NOTE**
These notes give you useful additional information on the installation or operation.

 These test strips for blood glucose tests correspond to the requirements of the EU guidelines 98/79 for in vitro diagnostic devices and bear the CE marking (conformance label) CE 0483.

 In vitro diagnostic medical device (for external use only)

 Biohazard  Contains sufficient for <n> tests

 Do not re-use  LOT number

 Use by  Catalogue number

 Storage temperature limitation  Manufacturer

 IMPORTANT NOTE

Read this leaflet and the instruction manual for the GlucoDock® blood glucose module carefully before using the test strips. If you still have any questions or require assistance, contact the MEDISANA service centre.

INTENDED USE

The GlucoDock® test strips are used with the GlucoDock® blood glucose module from MEDISANA to determine the blood glucose level in fresh capillary blood. They are intended for external use only (in vitro diagnostics). The GlucoDock® blood glucose module can be used at home, as well as professionally.

SAFETY INFORMATION ON STORAGE AND USE

- Do not store the test strips in the refrigerator or freezer.
- Store the test strips in a cool, dry place between 2 °C and 30 °C (35,6 °F - 86 °F).
- Only store the test strips in their original container. Never put them in a new or different container.
- Keep the test strips away from direct sunlight and any other heat source.
- Protect the test strips from high humidity.
- Write down the date you first opened the container on the test strip container.
- The test strips can no longer be used 3 months after opening. Dispose of the remaining test strips in the container together.
- Use each test strip directly after removing it from the container. Close the container immediately after removing a new test strip. Keep the lid sealed at all times.
- A test strip must be used within three minutes of removing it from the container.
- Only touch the test strip with dry, clean hands. Protect the test strip from getting contaminated.
- Never bend, cut or modify the test strips.
- Only use GlucoDock® test strips.

SAFETY INFORMATION FOR PATIENTS

- Only for use with in vitro diagnostics (for external use only). Do not put in the mouth or swallow.
- Never use the test strips after the expiry date.
- The test strips are only intended to be used once. Never use them twice.
- Medical staff and any others who use this system on more than one patient should be aware that all products or items that come into contact with human blood should be treated as if they are capable of spreading viral diseases, even after cleaning.
- Consult your doctor if you notice symptoms which do not correspond to your test result, despite having followed all the instructions in the manual for the module.
- You may only adapt the procedure for using products at home and self-monitoring, if you have first received the appropriate training to do so.

RESTRICTIONS

- GlucoDock® blood glucose test strips are for use with a fresh sample of capillary blood. Do not use serum or plasma.
- **Do not use the test strips for tests on newborn babies.**
- Extreme humidity can affect the test results. A relative humidity of more than 90% can lead to incorrect readings.
- The operating temperature of the system is between 10 °C and 40 °C (50 °F - 104 °F). If the temperature is not within this range, incorrect test results can occur.
- Haematocrit: the test results will not be affected by a haematocrit level of between 30% and 50%. A haematocrit level of below 30% can cause an incorrect (too high) test result, a haematocrit level of above 55% can cause an incorrect (too low) test result. If you do not know what your haematocrit level is, consult your doctor.
- Metabolites: the affect of lowering agents depends on the concentration. The following substances may affect the test result:
 - Azetaminophen $\geq 1,00$ mmol/L or 15,0 mg/dL
 - Gentisic acid $\geq 0,50$ mmol/L or 8,0 mg/dL
 - L-dopa $\geq 0,50$ mmol/L or 10,0 mg/dL
 - Dopamine $\geq 0,70$ mmol/L or 13,0 mg/dL
 - Methyldopa $\geq 0,12$ mmol/L or 2,5 mg/dL
 - Uric acid $> 0,40$ mmol/L or 14,0 mg/dL
- Patients undergoing oxygen therapy treatment may receive imprecise test results
 - The test strips can be used at altitudes of up to 3,048 m without having any impact on test results.
- Incorrect test results can occur if the patient is severely dehydrated, suffers from very low blood pressure or is in a state of shock, as well as in cases of hypoglycaemia or hyperosmolarity (with or without ketosis).
- Lipaemia effects: very high levels of blood fat in the blood sample could have a negative impact on certain methodologies. To rule this out, the patient should have his basic level checked in a clinical laboratory test by his doctor before performing the home test. Thereafter, the basic blood glucose level of the patient must be regularly checked and redetermined if necessary.

- Increased cholesterol levels and triglyceride can lead to incorrect test results in certain circumstances.
- Studies have shown that electromagnetism affects the electrical performance and precision of medical equipment and could therefore lead to an incorrect diagnosis.

MEASURING METHODS

The GlucoDock® blood glucose module was designed to quickly determine the blood glucose volume using electrochemical technology. This technology involves using disposable reagent strips. In a test, the electrical current is measured which is caused by the reaction of the glucose (glucose oxidase) with the reagent (Aspergillus niger) of the strip. The module measures this current and shows the corresponding blood glucose volume. The strength of the current caused by the reaction depends proportionally on the amount of glucose in the blood sample. In this way, the concentration of blood glucose is determined and is shown in the display of the GlucoDock® module within approx. 5 seconds.

CONTROL TEST TO ENSURE PRECISE RESULTS

A control test should be performed to ensure precise results of the GlucoDock® blood glucose module in the following cases:

- You are using the module for the first time.
- You have a new container of test strips.
- You have already repeated a test, and the results are still lower or higher than you expected.
- You suspect the module or test strips might not be working properly.
- The module has been dropped.
- The test strips have been exposed to a temperature which does not lie within the specified range for the storage conditions (2 °C - 30 °C or 35,6 °F - 86 °F).
- Perform a control test at least once a week.

A control test with control solution is performed to check the precision of the module and test strips. Proceed as described in the instruction manual of the GlucoDock® blood glucose module. The MEDISANA control solution is only to be used with the GlucoDock® blood glucose module from MEDISANA. Other control solutions may cause incorrect test results. If the result of a test with control solution lies within the range printed on the test strip container, this confirms that the module and the test strips are functioning perfectly. If the test result does not lie within the specified control range, repeat the control test. If the test result is still not within the permissible control range, check the standard procedure for the module, control solution and test strips according to the instruction manual. If incorrect results persist, contact the service centre.

Test results which deviate from the control range might be caused by the following:

- The expiry date of the control solution has been exceeded or the control solution is contaminated.
- An error occurred while performing the test.
- The module is not working correctly.
- The expiry date of the test strips has been exceeded or the test strip is damaged.

Do not use the module to measure your blood glucose level, until you have reached a control test result within the control range.

 **NOTE**
The control range on the test strip container can change with each new container. Always observe the current levels on the label when comparing.

TAKING A BLOOD SAMPLE

The GlucoDock® blood glucose module is intended for use with fresh capillary blood. A blood sample must be used in the test directly after drawing it. For a blood glucose test using the GlucoDock® module, a minimum blood sample of 0,6 µL is required.

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